

**Product Name :**  
Levelling Staves

**Product Code :**  
NLAB-SURVYLAB180005



**Description :**

Levelling Staves

**Technical Specification :**

**Specification:**

Level Staff is a graduated aluminum rod, used with a leveling instrument to determine the difference in height between points or heights of points above a datum surface. It cannot be used without a leveling instrument. a range of levelling staves including E-type, bar coded, standard invar and industrial Invar Levelling Staves. Most of the staves are sectional and so can be adjusted in length to allow for easy storage and transport. The aluminium staves adjust with telescopic section inside each other, whilst the wooden versions use sliding slip joints or connections.  
Sizes: 4/5/6 metres

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003