

Product Name : Lever Kit	Product Code : NLAB-EDUSCHOOLEQ174057
Description :	
Lever Kit	
Technical Specification :	
To illustrate principles of lever and moment. Five plywood beams 600 x 44 x 3 mm with ruled graduations at every 38 mm and a shallow groove across its centre for fulcrum point, with 4 wood fulcrums and forty 25 mm square metal blanks.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003