

Product Name : Light Fastness Tester	Product Code : NLAB-TESTINGLABM61002
--	--

--

Description :

Light Fastness Tester

Technical Specification :

Features:

Water Circulating Pump: Separate from Main Unit

Tubing: Flexible Transparent Type

Timer: Digital Timer with preset arrangement

Buzzer: Attached with timer

Water: Use of distilled water only in the water circulation unit.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003