


Product Name : Light Meter	Product Code : NLAB-TESTINGLABM19002
--------------------------------------	--



Description :

Light Meter

Technical Specification :

Technical Specifications

Measurement : 0 ~ 20,000 Lux (2 ranges)
Display : 13mm LCD, 3 ½ digits
Maximum indication : 1999
Resolution : 1 Lux and 10 Lux
Accuracy : ±(5% +2d)
Separate light sensor to measure light at an optimum position
Pocket size, easy to carry out & operation
Built in low battery indication
LSI circuit for high reliability
power : DC 9 volt battery

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003