

Product Name :
Light Sensor

Product Code :
NLAB-VOCATIONALTVET28012



Description :

Light Sensor

Technical Specification :

Resolution:

0-600 lux: 0.2 lux

0-6000 lux: 2 lux

0-150000 lux: 50 lux

We are well-known manufacturers, OEM suppliers of Light Sensor for Data Acquisition System. Contact us for high quality Light Sensor for Data Acquisition System for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com