

Product Name : Light Turbo 5 Watt Boost	Product Code : NLABMEDICAL-M32100017
Description : Light Turbo 5 Watt Boost India, Light Turbo 5 Watt Boost manufacturer India, suppliers, Medical Light Turbo 5 Watt Boost India, Hospital Light Turbo 5 Watt Boost India	
Technical Specification : Spectrum: 450-470nm Dimension: D28 X L222 mn Weight: 8 mm angled optic fiberg	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003