

**Product Name :**  
Linear IC Trainer

**Product Code :**  
NLAB-ELECTRICALAB210035



**Description :**

Linear IC Trainer

**Technical Specification :**

**SCOPE OF LEARNING:**

Study of Linear Circuits and Its applications.

**TECHNICAL SPECIFICATIONS:**

**Power Supplies:**

DC Power Supply IC Regulated +15V DC, 150mA (Fixed).

DC Power Supply IC Regulated 0-5V DC, 150mA (Variable).

DC Power Supply IC Regulated 0-12V DC, 150mA (Variable).

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

741 Operational Amplifiers: 2 Nos.

Resistors

Capacitors

Zener, Led and Diodes

Variable Resistor 10K

Control Switches: 2 Nos.

ZIF Socket 16Pin and 20Pin

Bread Board 1Nos.

#### SALIENT FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Connections are brought out through 2mm Brass Sockets.

Patch Cords 2mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

#### OPTIONAL ACCESSORIES:

Function Generator

CRO

Digital Multimeter

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003