

Product Name :
Linear Slide Staining Machine

Product Code :
NLAB-LABORATORYEM28013



Description :

Linear Slide Staining Machine

Technical Specification :

It has 14 stations, 1000ml capacity of single cylinder, with 1 to 59 minute(s) working time of single cylinder, with 72 glass slides, automatic cleaning and drying functions.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003