

**Product Name :**  
Linear Slide Staining Machine

**Product Code :**  
NLAB-LABORATORYEM28013



**Description :**

Linear Slide Staining Machine

**Technical Specification :**

It has 14 stations, 1000ml capacity of single cylinder, with 1 to 59 minute(s) working time of single cylinder, with 72 glass slides, automatic cleaning and drying functions.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003