

Product Name :
Linear Tissue Processor

Product Code :
NLAB-LABORATORYEM28014



Description :

Linear Tissue Processor

Technical Specification :

Battery Backup of upto 4 hrs.

Can be run both in manual and automatic mode.

Advance safety features to ensure the safety of tissue.

Program can be delayed start.

Easy to use with the help of touch screen.

Technical Specifications :-

Configuration: Linear

No. of stations: 12

Station capacity: 1000ml 2000ml

Paraffin stations: 3

Reagent stations: 9

Multi program memory: Yes

Temp. range (Paraffin): Ambient -80 deg C

Delay start: Yes

Battery backup: Yes

Safety device: Yes

Touch screen: Yes

Fume exhaust fan: Yes

Chamber illumination: Yes

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003