

Product Name :
Linen Change Trolley

Product Code :
NLABMEDICAL-M2620016



Description :

Linen Change Trolley

Technical Specification :

Specification :-

Frame made of round tube.

The Linen Change Trolley is divided into 2 parts.

One part has three shelves.

The other part has a good quality superior canvas bag.

The complete soiled linen trolley is mounted on 10 cms castors.

We are Manufacturer, Exporters and Suppliers of Linen Change Trolley at tender and export competitive prices from India.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments