

Product Name :
Liquid Limit Device (Motorized)

Product Code :
NLAB-TESTINGLABM109008



Description :

Liquid Limit Device (Motorized)

Technical Specification :

Same as above but fitted with a motor geared down to give approximately 120 rpm. Suitable for operation on 230 V A.C., Single Phase, 50 Cycles.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003