

**Product Name :**  
Liquid Limit Device (Motorized)

**Product Code :**  
NLAB-TESTINGLABM109008



**Description :**

Liquid Limit Device (Motorized)

**Technical Specification :**

Same as above but fitted with a motor geared down to give approximately 120 rpm. Suitable for operation on 230 V A.C., Single Phase, 50 Cycles.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003