

Product Name :
Liquid Limit Penetrometers

Product Code :
NLAB-SOILELAB140012



Description :

Liquid Limit Penetrometers

Technical Specification :

Penetration cones and sample cups has to be ordered separately.

Weight approx.: 8.5 kg

Penetrometer accessories:

Penetration test cone 30° angle

Cone test gauge to check the condition of the cone

Penetration sample cup 55 mm dia., 40 mm deep

Penetration sample cup 75 mm dia., 50 mm deep

Penetration test cone 60° angle, 60 g weight

Cone test gauge for cone 22-T0029/4

Penetration test cone 30° angle, 100 g weight

Penetration test cone 30° angle, 400 g weight

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003