

**Product Name :**

Liquid-Liquid Extraction Pilot Plant With Raschig Ring Packing

**Product Code :**

NLAB-ENGINEERINGLB36024

**Description :**

Liquid-Liquid Extraction Pilot Plant With Raschig Ring Packing

**Technical Specification :**

Feed tank with a capacity of 30 litres for the heavy phase material.

Feed tank with a capacity of 30 litres for the light phase material

Collecting tank with a capacity of 30 litres for the extracted product

Packed extraction column 600mm, with Raschig rings.

Electronic indicator for agitator speed.

Top separator.

Bottom separator.

Emergency stop button

2x metering pumps with a flow rate of 0 to 40 litres per hour.

Collecting tank with a capacity of 30 litres for the refined product

The electrical control panel in the main cabinet

We are well-known manufacturers, OEM suppliers of Liquid-Liquid Extraction Pilot Plant With Raschig Ring Packing for Chemical Engineering Lab Equipments. Contact us for high quality Liquid-Liquid Extraction Pilot Plant With Raschig Ring Packing for Chemical Engineering Lab Equipments for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)