

<b>Product Name :</b> Liver Infusion Powder	<b>Product Code :</b> NLAB-CHEMICALSIND20012
<b>Description :</b>	
Liver Infusion Powder	
<b>Technical Specification :</b>	
Liver Infusion Powder	
<p>Liver Infusion Powder is suitable for vaccine manufacturing and bacteriological culture media, for cultivation of fastidious organisms. Liver infusion powder is a dried infusion of fresh liver tissues of healthy animals. It is carefully digested with enzymes and dehydrated under controlled conditions to retain all nutritive values of liver.</p> <p>Solubility: Freely soluble in purified/distilled water, insoluble in alcohol.</p> <p>Storage: Store below 30°C. Use before expiry date on the label.</p> <p>Reaction: Reaction of 2% w/v aqueous solution at 25°C.</p> <p>Clarity: 1% w/v aqueous solution is clear without any haziness after autoclaving at 15 lbs pressure (121°C ) for 15 minutes.</p> <p>Appearance: Brownish yellow to brown Homogenous Free flowing powder ,having Characteristic odor but not putrescent.</p>	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003