

Product Name :

Load Frame for Universal Compression Test Of Rocks 12 Speed 200 Kn

Product Code :

NLAB-COMPSLAB280001

**Description :**

Load Frame for Universal Compression Test Of Rocks 12 Speed 200 Kn

Technical Specification :

The equipment consist of the following replaceable parts :

- Load frame, 200 kn capacity 12 speed from 0.0064 to 1.25 mm/min
- Proving ring, 200 kn with calibration certificate
- Dial gauge 25mm travel, 0.01mm least count.
- Platen set as per ASTM 2938 requirements

- Platen set as per ASTM 2938 requirements
- Horizontal clearance 300 mm
- Vertical clearance 750 mm
- Maximum platen diameter 198 mm
- Maximum platen travel 100 mm
- Specimen diameter 38 mm to 100 mm
- Can test AX, BX, NX rock specimen upto 100 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003