

Product Name :
Load Ring Penetrometer

Product Code :
NLAB-SOILELAB140022



Description :

Load Ring Penetrometer

Technical Specification :

Used for measuring the bearing strength and compaction degree of soils. The apparatus consists of a "T" shaped handle connected to a load ring 1 kN (100 kgf) cap., with max load pointer, and of an extension rod with five 100 mm graduations. The 30° end cone has an area of 645 mm² (1 sq. in). Supplied complete with calibration chart.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003