

**Product Name :**  
Locket Compass

**Product Code :**  
NLAB-ENGINEERINGLB87002



**Description :**

Locket Compass

**Technical Specification :**

Highly cherished among antique collection lovers.

Locket Compass is creatively designed by our innovative designers using optimum grade brass with the help of latest carving methodologies.

The compass provided by us can easily grab attention due to its attractive design and seamless finish.

Additionally our quality auditors check this Locket Compass on distinct parameters and can be availed from us at highly discounted prices within estimated time frame.

Features:-

Stunning look

Gleaming surface

Lively appearance

Eye-catching design

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003