

Product Name :
Logic Gate Trainer Using Cmos IC

Product Code :
NLAB-ELECTRONICSAB220004

Description :

Logic Gate Trainer Using Cmos IC

Technical Specification :

SCOPE OF LEARNING:

Logic Gates Using CMOS IC

TECHNICAL SPECIFICATIONS:

Power Supplies:

DC Power Supply +5V,150mA

Operated on Mains power 230V, 50Hz +10%

Technical specifications:

Eight independent logic level inputs to select High / Low TTL levels, each with a LED to indicate high / low status and termination.

Eight independent buffered logic level indicators for High / Low status

Circuit with OR, AND, Not, Nand, Nor, Ex-Nor & Ex-OR Gate given on board with Multi Operation

One BCD to 7 Segment Decoder given on Board.

One Bread Board is given on board.

SALIENT FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Instruction manual.

Connections are brought out through 2mm Colored Sockets.

Patch Cords 2mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003