

Product Name :
Logic Gates Circuit Trainer

Product Code :
NLAB-ENGINEERINGLB13001



Description :

Logic Gates Circuit Trainer

Technical Specification :

Logic Gates Circuit Trainer Lab Equipments Manufacturer, Suppliers & Exporters The system combines simple, easy to use, ICs for various gates and flip-flops, power supply, clock input and output state with a versatile solder less bread board area..

Student's new to digital electronics can implement logic circuits in a matter of minutes on the bread board area. Unit can be used by both the absolute beginner, and in advanced project work. The many outstanding features of the logic trainer, combined with its ease of use and robust housing, make it the first choice for those wishing to introduce students to digital electronics for the first time. Practical experience on this board carries great educative value for Science and Engineering Students

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com