

Product Name :
Logic Trainer

Product Code :
NLAB-TECHNICALAB17003



Description :

Logic Trainer

Technical Specification :

The logic trainer is designed for the logic beginners to enhance the comprehension of basic logical theory. The design is easy to operate and easy to understand. It is equipped with various kinds of basic logic gates, debounced logical switches, LED indicators, DC power supply with short circuit protection, pulse generator and solderless breadboard.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003