

**Product Name :**  
Los Angeles Abrasion Test Machine

**Product Code :**  
NLAB-TESTINGLABM55003

**Description :**

Los Angeles Abrasion Test Machine

**Technical Specification :**

Features

The Los Angeles Abrasion Test Machine Is Used As Quality Control Tool For Aggregates, Including Crushed Rock, Gravel And Slag

The Testing Method Determines The Degradation When Subjected To Abrasion, Attrition, Impact And Grinding

Heavy Steel Cylinder And Base Frame

Dust Proof Cover Plate

The Machine Comply With ASTM C131, EN1097-2, AASHTO T-96 International Standard

Specifications

Diameter Length : 710 X 510 Mm

Drum Rotation Speed : 30 ~ 33 Rpm

Counter : 5 Digits, Adjustable, Automatic Shut Off, Push Button Reset & Adjustment

Steel Balls : Ø48 Mm

Qty. Of Steel Balls : 12

Steel Tray For Specimen Unloading

Power Supply : 220V AC, 50/60Hz, Single Phase

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003