

Product Name :
Los Angeles Machine With Steel Bowl

Product Code :
NLAB-EDUCATIONTVET15038



Description :

Los Angeles Machine With Steel Bowl

Technical Specification :

internal diameter 711 mm x 508 mm wide. Rotational speed 31 to 33 rpm.

Standards: ASTM C131, CNR N ° 34, AASHTO T96, EN1097-2, UNE 83116, NLT 325, UNI 8520

Abrasive load set of 12 (according to EN and NF)

Pallet truck, capacity 2.5 T

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003