

**Product Name :**  
Losses in Pipes Friction Apparatus

**Product Code :**  
NLAB-SCIENTIFICIND13001



**Description :**

Losses in Pipes Friction Apparatus

**Technical Specification :**

Pipe Test Section: Dia  $\frac{1}{2}$ ", Length: 1m, Material G.I. Dia  $\frac{3}{4}$ ", Length: 1m, Material G.I.

Water Circulation: FHP Pump, Crompton make.

Flow Measurement: Using Measuring Tank with Piezometer Capacity 25 liters.

Sump Tank: Capacity 50 Liters.

Stop Watch: Electronic

Control Panel Comprises of: Standard make On/Off Switch, Mains Indicator, etc.

Tanks will be made of Stainless Steel.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments