

Product Name : M8B Multipara Monitor	Product Code : NLABMEDICAL-M3360005
--	---

--

Description :

M8B Multipara Monitor

Technical Specification :

M8B Multipara Monitor

Lightweight, compact, portable.

10 inch wide color TFT LCD.

Nellore Oximax/EDAN SpO2.

SpO2 pulse-tone modulation (Pitch Tone).

Multi-leads ECG waveforms (7 leads) display.

Arrhythmia analysis, multi-lead S-T segment.

Built-in rechargeable Lithium-ion battery.

Network capability.

Powerful data storage capacity (96 hours graphic and tabular trends of all parameters 500 NIBP measurement, 60 alarm events).

Optional accessories for adult, pediatrics and neonate.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003