

Product Name :
Mechanical Double Action deep Drawing Presses

Product Code :
NLAB-MACHINELABWK14003



Description :

Mechanical Double Action deep Drawing Presses

Technical Specification :

Mechanical Double Action deep Drawing Presses .

PRECISION

BEARING

487.5 x 508 x 830
1000
650

mm mm

PRE~~SSURE~~
IN TONS

600~~ANK~~

ST~~RUCKE~~

Mm~~Min~~

H0~~5725~~"

48~~17365~~ / 600

mm mm

12~~71AX~~

0~~5038E~~

mm mm

17~~2'AP~~

C1~~300CH~~

M~~11DIA~~

F~~106~~HT

2000

WITHOUT

MOTOR

**M35HP
HPM**

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003