

Product Name :
Machine for Accelerated Polishing

Product Code :
NLAB-EDUCATIONTVET15041



Description :

Machine for Accelerated Polishing

Technical Specification :

Description:

Used to assess the resistance of road aggregates to polishing under the effect of traffic.

The machine includes a test tube wheel, capable of receiving 14 test tubes, which is rotated and on which a rubber tire is applied.

The wheels are stacked and the abrasive sand is poured between them

Characteristics:

Possibility to place 14 samples

Power supply: 220 V - 50 Hz

Weight: 175 kg

Accessories:

Spare mold

Mold cover

Emery grain (coarse), 6 kg

Emery powder (fine), 6 kg

Test stones, 25 kg

Reference stones (Criggion), 25 kg

7.2x40 mm flattening screen

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003