

Product Name : Magnetic Needle	Product Code : NLAB-EDUSCHOOLEQ1730207
Description : Magnetic Needle	
Technical Specification : A 62mm long, magnetic needle made from carbon steel with a brass cup bearing for low friction pivoting when used with the needle stand. Features:- Ideal for the study of compasses and magnetic fields. Accurately balanced. Length 62mm.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003