

Product Name : Magnetic Needle	Product Code : NLAB-EDUSCHOOLEQ1730207
Description :	
Magnetic Needle	
Technical Specification :	
A 62mm long, magnetic needle made from carbon steel with a brass cup bearing for low friction pivoting when used with the needle stand.	
Features:-	
Ideal for the study of compasses and magnetic fields.	
Accurately balanced.	
Length 62mm.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003