

Product Name :
Magnetic Stirrer cum Hot Plate

Product Code :
NLAB-LABORATORYEM42011



Description :

Magnetic Stirrer cum Hot Plate

Technical Specification :

Magnetic Stirrer cum Hot Plate

Features:-

Body made of CRC sheet
duly powder coated with top plate made of Stainless Steel.

A heavy duty PMDC motor
imparts a higher initial torque for smooth and vibration free start at lower
speeds as well as at higher speeds.

Variable Speed of up-to
1200 RPM is controlled through a step less speed controller.

For heated magnetic
stirrers the electrical/heat energy is controller through a simmer stat.

To work on 220/230 Volts
AC Supply.

Technical:-

Capacity Plate Dimensions (Plate) Ø135	
2.06" ltrs.	
6.08" ltrs.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003