

Product Name :
Magnetizing and Demagnetizing Solenoid

Product Code :
NLAB-EDUSCHOOLEQ174076



Description :

Magnetizing and Demagnetizing Solenoid

Technical Specification :

A helical coil 310 mm long of DCC copper wire mounted on a wooden board, number of turns marked, with one meter connecting cable, fuses, push switch compass & soft iron rod 350 x 10 mm. Internal dia of coil 25 mm approx. For use on 230 V mains supply.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003