

**Product Name :**  
Manual Batch Distillation Pilot Plant with Data

**Product Code :**  
NLAB-ENGINEERINGLB36001



**Description :**

Manual Batch Distillation Pilot Plant with Data

**Technical Specification :**

Manual Batch Distillation Pilot Plant with Data Lab Equipments Manufacturer, Suppliers & Exporters

Manual Batch Distillation Pilot Plant With Data Acquisition

- Distillation column height 600 mm with 7 bubble cap trays.
- Reboiler with 4 litres capacity complete with 1600 W electric heater
- Condenser with reflux head
- Tank for collection of distillate with 1 litre capacity
- Tube-in-tube heat exchangers manufactured from AISI 304 stainless steel

- 6 thermocouples with probes protected by stainless steel shield sleeves
- 6 electronic temperature indicators
- Thyristor unit for controlling the power to the reboiler electric heater.
- 2 programmable timers for the control of reflux ratio
- Variable-area flow meter made of stainless steel with a range 20 to 250 litres / hour and 4-20mA output signal.
- Vacuum gauge range -1 to 0 bar
- Electrical control panel in main cabinet
- PC interface contained in main electrical control panel
- Data acquisition software for Windows
- Vacuum circuit with trap
- Pipe work and valves manufactured from stainless steel
- Emergency stop button

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003