

Product Name :
Manual Bioethanol Production Plant

Product Code :
NLAB-ENGINEERINGLB36034



Description :

Manual Bioethanol Production Plant

Technical Specification :

Insulation of Rockwool and sheet steel

Thermocouple with Pt100 probe

Sensor for measuring pH

Safety valve

A distillation column with bubble cap trays and two thermocouples with PT100 probes

Multifunction reactor with 13 litres capacity complete with:

Condenser manufactured from stainless steel and two thermocouples with Pt100 probes.

Electrical control panel in main cabinet.

Emergency stop button.

Stirrer with variable speed and mechanical seal

Baffles

Sight glass

Cooling and heating coil

We are well-known manufacturers, OEM suppliers of Manual Bioethanol Production Plant for Chemical Engineering Lab Equipments. Contact us for high quality Manual Bioethanol Production Plant for Chemical Engineering Lab Equipments for schools lab, college lab, universities, research labs, various teaching and

workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003