

Product Name :
Manual CBR Press

Product Code :
NLAB-EDUCATIONTVET13024



Description :

Manual CBR Press

Technical Specification :

Capacity 50 kN

Supplied with 50 kN torque ring and comparator

Penetration speed display

Torque ring 30 kN (corresponding to a range of average CBR value)

Torque ring 2 kN (CBR value <8%)

Torque ring 10 kN (CBR 8 to 40%)

Device for maintaining the maximum value

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com