

Product Name :
Manual Cooling Tower

Product Code :
NLAB-ENGINEERINGLB36041



Description :

Manual Cooling Tower

Technical Specification :

Water pump with a maximum flow rate of 3 m³per hour and head of 5 m H₂O

6x digital thermometers with LCD displays

Flow meter with a range of 20 to 200 liters per hour.

Rectangular tower manufactured of Plexiglas with a height of 550 mm.

Centrifugal fan with a maximum flow rate of 1,340 m³per hour and head of 80 mm H₂O.

Tank manufactured from stainless steel for hot water with a capacity of 9 litres complete with 3 x 500 watts electric heaters.

The electrical control panel in the main cabinet

Emergency stop button

The 3x different packaging of 92,131 and 235 m²/m³ easily interchangeable in the tower

Additional water tank manufactured from Plexiglas with a capacity of 1 litre

We are well-known manufacturers, OEM suppliers of Manual Cooling Tower for Chemical Engineering Lab Equipments. Contact us for high quality Manual Cooling Tower for Chemical Engineering Lab Equipments for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003