

**Product Name :**  
Manual Cryostat Microtome

**Product Code :**  
NLAB-LABORATORYEM28015



**Description :**

Manual Cryostat Microtome

**Technical Specification :**

Easy to operate control panel.

Digital display of set and actual temperature.

Electronic temperature control.

Spacious, Stainless steel cryo chamber i.e. Easy to clean and disinfect.

Anti glare illumination.

Heated glass window for clear visibility.

Automatic defrosting system.

Full body stainless steel construction.

Built-in drainage system.

Easy mobility with heavy duty castors.

Hermetically sealed and thermostatically control Refrigeration.

Technical Specification :-

Operation: Manual

Chamber temp. range: 0 to -30°C / -40°C

Section thickness range: 1-20 microns

Trimming thickness range: NA

Specimen size (section): 40 x 55mm

Horizontal specimen feed: 22mm

Vertical specimen stroke: 55mm

Trimming: Manual

Programmable start / stop: Optional

Cooling down time: 60 mins.

Freeze shelf temperature: NA

Pettier unit of F.S: NA

Automatic defrost: Optional

Motorized sectioning: NA

Knife holder adjustments: Adjustable

Knife: 120mm

Refrigeration: Single compressor

Chamber size: 241x432x610mm

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003