

Product Name :
Manual Hand Shear

Product Code :
WeldingLab0025



Description :

Manual Hand Shear

Technical Specification :

12" Manual Hand Shear.

Blade length: 12"

Hardened high carbon steel blade

Heavy Duty steel frame

Maximum Cut Flat Steel: 2-3/4 (in) x 1/4 (in) (70 x 6mm)

Maximum Cut Rod Steel: 1/2 (in) (13mm)

Maximum Cut Thickness: 1/4 (in) (6mm)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773