

Product Name :
Manual Lensometer

Product Code :
NLABMEDICAL-M3230011

Description :

Manual Lensometer India, Manual Lensometer manufacturer India, suppliers, Medical Manual Lensometer India, Hospital Manual Lensometer India

Technical Specification :

Cordless operation and less power consumption

Compact, light weight and handy instrument

Simple and easy, scale reading system

Specifications:

Type: Internal/External Reading

Target: Cross & Corona., 360 Degree.

Vertex Power Scale: 0 to +_ 24 diopters

Prism Diopter Scale: 0 ^ to 5(1 ^step)

Cylindrical Axis: 0 Degree to 180 Degree (5 Degree step)

Lens Diameter: Upto 90 mm

Eyepiece Adjustment: +3 to -5 diopters, 0 to 90 degree, freely adjustable

Light Source: Two 1.5 V AA batteries (Also, Usable with 3v AC adapter)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

