

Product Name :
Manual Rotary Microtome

Product Code :
NLAB-LABORATORYEM28017



Description :

Manual Rotary Microtome

Technical Specification :

The Manual Rotary Microtome includes a knife holder capable of handling both high and low profile blades. The knife holder has smooth lateral movement to utilize the full length of the cutting blade and slides easily onto the dovetail base. For added safety, the knife holder comes equipped with an attached knife guard to cover the blade when the instrument is unattended. The knife guard also contains an ejecting pin mechanism for safe and simple removal and disposal. The knife clearance angle, 0 – 10°, can be adjusted on the knife holder base. A waste tray slides effortlessly under the front of the unit and is easily removed for cleaning.

Features:

Universal cassette block holder

with self-contained level

Versatile lock on hand wheel

Full lateral movement of knife holder

Safety knife guard/blade ejector

Easy-to-clean waste tray

Compact, sleek and sturdy design

Specifications:

Feed advance: 28mm

Vertical stroke: 60mm

Thickness 0.5 – 60µm

Dimensions: 17 ¼" x 21 ¼" x 11 ½" (440 x 540 x 290mm)

Weight: 51.4 lbs (23.3kg)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003