

**Product Name :**  
Manual Stirred Continuous Reaction Pilot Plant With Reactors In Series

**Product Code :**  
NLAB-ENGINEERINGLB36031



**Description :**

Manual Stirred Continuous Reaction Pilot Plant With Reactors In Series

**Technical Specification :**

Collecting tank with a capacity of 20 litres

Tank with a capacity of 2 litres for storing the neutralizing solution

Metering pump for the neutralizing solution with a maximum capacity of 1.5 litres per hour

4x thermocouples with pt100 probes protected by stainless steel shield sleeves

4x electronic temperature indicators

2x electronic conductivity meters

2x reactors with a capacity of 3 litres

2x feeding tanks with a capacity of 10 litres

The electrical control panel in the main cabinet

Emergency stop button

2x metering pumps with a flow rate of 0 to 8.4 litres per hour

2x stirrers with variable speed for the reactors

2x electric heaters for the reactors

We are well-known manufacturers, OEM suppliers of Manual Stirred Continuous Reaction Pilot Plant With Reactors In Series for Chemical Engineering Lab Equipments. Contact us for high quality Manual Stirred Continuous Reaction Pilot Plant With Reactors In Series for Chemical Engineering Lab Equipments for schools

lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003