

Product Name :
Manual Tray Dryer Process Unit India

Product Code :
NLAB-ENGINEERINGLB36042



Description :

Manual Tray Dryer Process Unit

Technical Specification :

4x aluminum drying trays 400 mm x 300 mm.

Distillation column height 600 mm with 7 sieve trays.

Horizontal air duct with 450 mm x 450 mm window.

Axial fan with digital speed control 0 to 4 meters per second.

Electrical control panel in main cabinet.

Electrical heating element array 2,700 watts.

2x digital thermohygrometers with digital displays.

Emergency stop button.

Digital anemometer with digital display.

Electronic balance with range 0 to 10 Kgs and 0.1 g resolution.

Digital chronometer.

We are well-known manufacturers, OEM suppliers of Manual Tray Dryer Process Unit for Chemical Engineering Lab Equipments. Contact us for high quality Manual Tray Dryer Process Unit for Chemical Engineering Lab Equipments for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com