

Product Name : Manual Vacuum Pump With Watch	Product Code : NLAB-EDUSCHOOLEQ1730169
Description :	
Manual Vacuum Pump With Watch	
Technical Specification :	
Made of plastic, extremely light-weight and user-friendly because of its ergonomic design.	
Ideal device for generating vacuum independent of mains supply.	
The pump reaches a difference. Against the external pressure of approx. Inflator	
Maintenance-free, self-lubrication and corrosion-resistant.	
Size: 23 x 18 x 4.5cm	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003