

Product Name :
Masses Hangers Iron

Product Code :
NLAB-EDUSCHOOLEQ174080



Description :

Masses Hangers Iron

Technical Specification :

Made of steel, chrome plated, designed to be used with various slotted masses..

Size: 100 g, 50 g, 20 g, 10 g

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003