

**Product Name :**  
Masses Hangers Steel

**Product Code :**  
NLAB-EDUSCHOOLEQ174081



**Description :**

Masses Hangers Steel

**Technical Specification :**

Made of iron, designed to be used with various slotted masses.

Size: 100 g, 50 g, 20 g, 10 g.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003