

Product Name : Material Storage Cabinet	Product Code : NLABMEDICAL-M3270008
Description :	
Material Storage Cabinet	
Technical Specification :	
Large space	
Safe storage	
Long durability	
Working Size: 3'x2'x5' 4'x2'x5'	
Size of HEPA Filters: 3'x1½'x6" 4'x1½'x6"	
No. of HEPA Filters: 1 1	
No. of Pre Filter: 1 2	
No. of Pre Filter: 3 3	
Whole Cabinet made of commercial board & laminated with sun mica.	
Whole Cabinet made of cold rolled steel sheet duly powder coated inside & outside for longer life durability.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003