

Product Name :
Maximum and Minimum Thermometers

Product Code :
NLAB-LBEMALAB230025



Description :

Maximum and Minimum Thermometers

Technical Specification :

A maximum–minimum thermometer. The scales are Fahrenheit on the inside of the U and Celsius on the outside. The current temperature is 23 degrees Celsius, the maximum recorded is 25, and the minimum is 15, both read from the base of the small markers in each arm of the U tube.

Ranges Available: -30° C to +50 ° C

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003