

**Product Name :**  
Maximum and Minimum Thermometers

**Product Code :**  
NLAB-LBEMALAB230025



**Description :**

Maximum and Minimum Thermometers

**Technical Specification :**

A maximum–minimum thermometer. The scales are Fahrenheit on the inside of the U and Celsius on the outside. The current temperature is 23 degrees Celsius, the maximum recorded is 25, and the minimum is 15, both read from the base of the small markers in each arm of the U tube.

Ranges Available: -30° C to +50 ° C

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003