

Product Name :
Maxwell Needle

Product Code :
NLAB-EDUSCHOOLEQ174089



Description :

Maxwell Needle

Technical Specification :

To determine the modulus of rigidity of metal wires or moment of inertia of rods. A hollow metal cylinder with pointer and chuck, two each solid and hollow metal cylinders, each of 50 mm length, with a metal wall bracket with a chuck.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003