

Product Name :
Meat Mincer

Product Code :
NLAB-VOCATIONALTVET27076

Description :

Meat Mincer

Technical Specification :

Capacity: 100kg/h

Voltage(V): 220 V

Weight(kg): 25kg

Material: Stainless Steel

Frequency: 50 Hz

Features:

Minimum maintenance

Excellent functionality

As per the industry standards

Longer working life

Seamless finish

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003